

Customer Supplied Hardware Requirements

Smartphone Requirements

mVisum, Inc. does not provide cellular hardware or mobile communication plans. It is the responsibility of the user to obtain appropriate equipment and mobile communication plans.

Minimum handheld smartphone requirements are as follows:

Blackberry Curve or 8800 series or superior with equivalent or superior screen size, pixel density, black depth, and with trackball navigation capability; minimum 1GB memory expansion card

iPhone 3G with operating system version 2.2.1 or later. mVisum will NOT be responsible for any device lock-outs should mVisum be attempted to be installed on unlocked/unauthorized (“jailbroken”) iPhone’s.

Smartphones that operate on the **Android** operating system (Nexus One, Motorola Droid, etc.) running Android v 1.6 or higher, with a minimum screen resolution of 360x480, 1GB or larger memory expansion card.

Windows Mobile Smartphone 5.0 or higher with touch screen capability, 128MB SD RAM/64MB Flash ROM, memory expansion slot with compatible 500MB or larger memory expansion card, screen size of at least 240x240 pixels with 16 bit color, GPRS data connectivity or superior, 150MHz processor or superior.

NOTE: Mobile phone plans should include unlimited data plan and SMS messaging plan commensurate with expected messages per month, 2.5 or 3rd Generation data connectivity is recommended.

NOTE: All mobile phone plans should include sufficiently large data plans and commensurate with expected messages per month. 2.5 G, 3G or 4G data connectivity is recommended for use over the cell phone network. Devices should be wi-fi capable for in-hospital use.

PC Requirements for mVisum Consoles and Physician Review Stations

mVisum, Inc. does not supply or sell computer hardware. It is the responsibility of the user to obtain appropriate computer systems or utilize existing computer systems and provide necessary network access.

Minimum system requirements for mVisum Consoles or Physician Review Stations are as follows:

Windows XP or later, 512 MB or higher RAM, 1GHz or higher processor, CD-ROM drive, Broadband access to the internet. Any installed firewall (within the computer or as part of the network) should allow access by the mVisum console software to mVisum servers.

Minimum Server Requirements

mVisum, Inc. does not supply or sell computer hardware. It is the responsibility of the facility to obtain appropriate computer systems or utilize existing computer systems and provide necessary network access.

Windows 2003 Standard Edition or Higher with QuadCore Xeon Processor or equivalent, 4GB or higher RAM, 250GB or larger hard drive space with RAID, SQL server 2005 Standard Edition or Higher with necessary connectivity to external cell phone data networks and SMS gateways, SQL Reports, IIS, and .Net 2.0.

mVisum can be deployed in a virtual server environment. Other software components may be needed depending upon integration requirements into hospital information systems.